

# iOptima<sup>INT</sup> Buying Guide




	1700704-IPAD Single Motor	1700705-IPAD Dual Motor	1700706-IPAD All-in-One
Internal Installation	✓	✓	✓
Restorative	✓	✓	✓
Rotary Endodontics WaveOne® Gold & RECIPROC® Ready	✓	✓	✓
Surgery/Implant Placement	✗	✗	✓

## SPECIAL BRACKETS

Each iOptima<sup>INT</sup> motor system includes a left/right bracket for installation onto most delivery units. For full integration into the following delivery units, order their corresponding brackets.

- ADEC 300.....REF 1307830-001 & 1304437-HRD
- ADEC 500.....REF 1307830-001
- ADEC 532B/533B (LEFT SIDE MOUNT ONLY).....REF 1308204-001
- ADEC 541 REAR DELIVERY UNIT .....REF 1900346-001

If you have a legacy delivery unit or are unsure, contact your local Bien-Air sales representative for additional custom bracket information.



Shown: REF 1307830-001 for ADEC 300 & 500 Units

S W I S S  M A D E

# iOptima<sup>INT</sup> Buying Guide



## ADD-ON KITS & ACCESSORIES

À la carte expansions and additions to existing iOptima<sup>INT</sup> systems or add a third attachment to a dual motor system. Each system can have up to 3 motors, or 2 motors with 1 scaler.

**MX2 PLUS LED KIT**  
REF 1700704-2MT

Additional motor for restorative & endo



**MX-i PLUS LED KIT**  
REF 1700704-MXI

Additional motor for surgery & implantology



**STERILE IRRIGATION PUMP**  
REF 1900200-001

Irrigation pump for surgical attachments



**PIEZO SCALER KIT (w/ LIGHT)**  
REF 1900300-IOPTIMA

**PIEZO SCALER KIT (w/o LIGHT)**  
REF 1900301-IOPTIMA



**iOPTIMA<sup>INT</sup> PUMP/SCALER RELAY BOARD**  
REF 1800026-001

For full connectivity and function of the Sterile Irrigation Pump and Piezo Scaler in the same system. Provided for FREE when both are purchased together.

Connect with us and schedule a FREE in-office demo!



<https://www.bienairusa.com/connect-with-us>

S W I S S  M A D E